

The Quick Discrimination of Nylon by ionRocket-DART®-MS

[BACKGROUND]

Nylon, the synthetic polymer and aliphatic polyamide, is one of the most widely used industrial polymers. Nylon is widely diffused due to its high rigidity and elasticity. nylon-6,10 and nylon-6,12 reflect similar tendencies in FT-IR spectrum. In addition, the melting temperatures of both polymers measured by DSC are approximately 220°C. This makes it difficult to discriminate between the two analogs.

[SAMPLES] nylon-6,10 and nylon-6,12 (e.g. toothbrush).

[METHODS]

The analyzing tool was the ionRocket heating system connected to the DART®-MS (Direct Analysis in Real Time - Mass Spectrometry). Small quantities of the samples were put on the POT and analyzed. The temperature was increased by 100°C per minute, from 30°C to 600°C.

[RESULTS]

TIC is shown in Fig 1. MS spectra measured at 300°C is shown in Fig 2.

In Fig 2, the monomers and dimers were mainly observed.

The comparisons of ionRocket-DART®-MS spectra made it possible to distinguish between nylon-6,12 and nylon-6,10.





Fig. 1 TIC of each sample.
R.T. 100°C/min 600°C




Fig.2 MS spectrum measured at 300°C
The preset temperature of DART®-SVP was 400°C
Ionization was DART® positive.

[TARGET] Material development, chemical industry, foreign material analysis