
Evaluation of degradation degree of lubricating oil by thermal desorption and pyrolysis combined with DART-MS (TDP/DART-MS) Chikako Takei, Kenichi Yoshizawa ; BioChromato, Inc.

Introduction

Lubricating oil is composed of base oil and additives. In order to analyze the base oils and additives of lubricating oils complicated pretreatment which takes a lot of time and effort were required generally. (Fig. 1). However for R&D, QC and market research, it is important to obtain the information on base oils and additives.

Recently, thermal desorption and pyrolysis/direct analysis in real time (TDP/DART)-MS^{1~3)} and Kendrick Mass Defect (KMD) analysis⁴⁾ is using for polymer analysis, respectively.



The purpose of this work :

Analysis directly and evaluate the degradation degree of lubricating oils.


Materials and Methods

Materials Automotive engine oil *Mileage : 0km, 1000km, 5000km

 Analytical methods TDP/DART-MS (Fig.3)

1µL of sample were put into the POT. Mass spectra were measured as the samples were heated.


Kendrick Mass Defect(KMD) analysis was used a "Spectra Scope (BioChromato)" software.




Results and discussions


In order to bird's-eye view of the whole mass spectra, heat map (horizontal axis: m/z, vertical axis: temperature) of each samples were shown in Fig. 4.

were detected.





Pyrolysis region ~ Base oil ~



Compounds with locally high strength at thermal desorption region (300 °C) and with repeating structure which seems to be polymer at pyrolysis region (450 °C)

Thermal desorption region ~ Additives ~

• A phenol type antioxidant, an amine type antioxidant, a salicylic acid type detergent dispersant were detected and determined from all samples(Fig. 5, 6). No significant difference in the amount ratio of additive components at 0 km and 1000 km running, but at 5000 km, salicylic acid type detergent and phenolic antioxidant was remarkably decreased(Fig. 5).

Conclusion

A combination of TDP/DART-MS and KMD analysis enables analysis both additives and base oil without any pretreatment.

Therefore, the combination can contribute to elucidate of the degradation mechanism, failure analysis, R&D, and quality control in the field of automobiles.

References

- C.Takei, K.Kinoshita, T.Nishiguchi, h.Shimada, K.Maeno, Y.Shida; 63rd ASMS Annual Conference, Poster ThP74 (2015)
- 2. R.B.Cody, C.Takei, H.Shimada, Y.Shida, A.Kusai; 64th ASMS Annual Conference, Poster MP545(2016)
- C.Takei, K.Yoshizawa, H.Ohtani; 65th ASMS Annual Conference, Poster WP605(2017)
- H.Sato, S.Nakamura, K.Teramoto, T.Sato; JASMS, 2014, 25(8), 1346-1355
- Minoru Fujita "Petroleum analysis", Petroleum Analysis Chemistry Lab (1992)

