

Analysis of Poly(butyl acrylate) by Py-GC/MS, MALDI-MS and Thermal Desorption and Pyrolysis combined with DART-MS

Chikako Takei ¹, Kenichi Yoshizawa ¹, Toshiji Kudo ², Hajime Ohtani ³; ¹ BioChromato, Inc., ² Bruker Japan K.K., ³ Nagoya Institute of Technology



Introduction

Pyrolysis-gas chromatography/mass spectrometry (Py-GC/MS) and matrix-assisted laser desorption/ionization (MALDI)-MS are widely used for polymer characterization.

Recently, thermal desorption and pyrolysis/direct analysis in real time (TDP/DART)-MS is using for polymer analysis¹⁻³.


Analytical purposes of polymer products range from raw material's lot survey to failure analysis. Thus, we are aiming at further development of polymer analysis technology, especially for base polymer constituting a polymer products. In this poster, we described our progressing using homopolymer as our first step.

The purpose of this work :

In the view of **repeating unit**, **molecular weight distribution** and **terminal functional groups**, clarifying the characteristics of each analytical methods by measurement data and general theory.

Experimental

◆ Material : Poly(butyl acrylate), PBA
(Sigma Aldrich)
*Mw ~99,000 (GPC)



◆ Analytical methods

(a) Py-GC/MS

Pyrolyzer: 600 °C
Carrier: He, 100 ml/min
Injection: Split 100 : 1 (1 ml/min at the column), 320 °C
Column: Frontier Lab UA-5 (30 m x 0.25 mm id, 0.25 μm)
Oven: 40 °C (2 min) - (20 °C/min) - 320 °C (13 min), Total 29 min
Detector: EI ionization
Quadrupole, scan range m/z 29-600
Sample: ca. 0.2mg

(b) MALDI-MS

Matrix: DCTB
Ionization accelerator: NaTFA
Solvent: THF
Laser: Smartbeam II (355 nm)
Analysis mode: Positive Reflector (accel. 19 kV, reflect. 21 kV)
Positive Linear (accel. 19.5 kV)

(c) TDP/DART-MS


TDP device: RT (1 min) - (100 °C/min) - 600 °C, Total 7min
Detector: DART Ionization (Ionization gas: He, 400 °C)
Q-TOF, Mass range m/z 50 - 2000
Sample: ca. 0.2mg

Kendrick Mass Defect (KMD) analysis was used a "Spectra Scope (BioChromato)" software.


Results and Discussions

Repeating unit


Py-GC/MS enables identification of polymer type using the available spectral library of pyrolysis products (ex. F-Search, Frontier Lab), (Fig. 1).



MALDI-MS enables identification of polymer type using the observed mass interval corresponding to the mass of the repeating unit or the elemental compositions calculated from the mass (Fig. 2).



TDP/DART-MS enables identification of polymer type without spectral library using the observed elemental compositions of the repeating unit (Fig. 3).



Conclusion

Regarding the measurement samples, Py-GC/MS and TDP/DART-MS enables measurement polymer materials such as vulcanized rubbers and crosslinked resins, but they disables grasp information about molecular weight distribution as-is. On the other hand, MALDI-MS enables measurement absolute molecular weight distribution as-is only if samples were soluble and have narrow distribution ($M_w/M_n < 1.1$). Additionally it enables analysis compositional characterization of polymers only if the spectra was measured (excluding vulcanized / crosslinked polymer, fluorine polymer and olefin polymer).


TDP/DART-MS enables individual measurement of low molecular weight pyrolysis products and oligomeric thermal desorption products, since it utilize a DART soft ionization technique combined with TDP device. Therefore, it give us valuable clue to elucidate chemical structure of unknown compounds (ex. Polymers which has wide distribution and compounds not listed in spectral library such as reacted and degraded products). Thus, it can contribute to further development of polymer analysis technology.

Our next step, we will clarify the characteristics of these analytical methods using copolymer, crosslinked polymer and degradation products.


Molecular weight distribution

Py-GC/MS disables observation of molecular weight distribution since their analysis target is pyrolysis products.

MALDI-MS enables observation of absolute molecular weight distribution, in case the sample has quite narrow distribution. This PBA has too wide to discuss about average molecular weight (Fig. 4, 5). In that case, SEC/MALDI-MS is effective (Fig. 6).




TDP/DART-MS enables measurement of molecular weight distribution of thermal desorption compounds (Fig. 7).




Terminal functional groups

Py-GC/MS disables identification of terminal functional groups of unknown polymer samples. But for a known terminal functional groups, quantitative analysis is possible.

MALDI-MS enables analysis how many terminal functional groups. KMD analysis, especially Remainder (R)KM plot⁴ is helpful, since the number of terminal functional groups are visualized clearly (Fig. 8).



TDP/DART-MS enables not only analysis how many terminal functional groups but also elucidation chemical structure of terminal functional groups using accurate mass and MS/MS spectra of thermal desorption products (Fig. 9). KMD and RKM plot⁴ are also helpful (Fig. 10).



References

- C. Takei, K. Kinoshita, T. Nishiguchi, H. Shimada, K. Maeno, Y. Shida; 63rd ASMS Annual Conference, Poster ThP74 (2015)
- R.B. Cody, C. Takei, H. Shimada, Y. Shida, A. Kusai; 64th ASMS Annual Conference, Poster MP645(2016)
- C. Takei, K. Yoshizawa, H. Ohtani; 65th ASMS Annual Conference, Poster WP605(2017)
- S. Tsuge, H. Ohtani, and C. Watanabe; Pyrolysis-GC/MS Data Book of Synthetic Polymers, Elsevier (2011) p. 92-93
- H. Sato, S. Nakamura, K. Teramoto, T. Sato; JASMS, 2014, 25(8), 1346-1355

Acknowledgment

I would like to express my very great appreciation to Dr. Hiroaki Sato(AIST) and Ms. Sayaka Nakamura (AIST) for their assistance with the collection of SEC/MALDI-MS data.