

## RAPID Slit Seal Instruction

Thank you for purchasing our RAPID Slit Seal.  
Please read this manual carefully before use and  
save this manual for your future use.

Please note that product specifications and packings  
may change without prior notice.

### 【Purpose of use】

This product can be used for research purposes only.

### 【Notes for use】

Handle this product carefully when unpacking from the bag.

Do not process this product in autoclave treatment.

### 【Storage and disposal of product】

Store at room temperature and in a place not exposed to direct  
sunlight. To dispose the product after use, comply with local laws  
and ordinances, and process properly according to the regulations  
of your facility.

## Instructions

\*Please wear gloves before use.


1

The film with a red indication is a release liner.


赤線の付いたシートが、剥離フィルムです

Peel off the release liner to expose the adhesive surface.

粘着面が露出します



2



Remove the release liner with a red indication and expose the adhesive  
surface.

Paste the RAPID Slit Seal on the 96-well plate and make sure  
holes of RAPID Slit Seal are correctly positioned on each well  
spot.

3



Remove the protective film before sampling or dispensing. If used  
without removing the protective film, the needle or the tip can not  
puncture the film. If the protective film stays on and is left at 100°C,  
the edge of the protective film may be curved.

## Manufacturer



BioChromato, Inc.

1-12-19 Honcho, Fujisawa, Kanagawa-ken 251-0053 JAPAN

TEL : +81-466-23-8382 FAX : +81-466-23-8279

E-mail : [infobc@bier.co.jp](mailto:infobc@bier.co.jp)

Website: <https://biochromato.com/>

Catalog Code	RSS-S96-10-S R80.120.10-S
	RSS-S96-50-S R80.120.50-S
Product name	RAPID Slit Seal Sterile 10sheets RAPID Slit Seal Sterile 50sheets
Bonding method	Pressure sensitive adhesive
Materials	PET, Silicone
Operating temperature limit	-80°C~100°C *
Size	W 80 mm × L 122 mm
Quantity	RSS-S96-10-S R80.120.10-S : 10 sheets/Box
	RSS-S96-50-S R80.120.50-S : 10sheets/Box × 5Boxes

\* The upper limit of the operating temperature is specified as 100 °C, and this is clarified by the evaluation test carried at  
100 °C for 30 minutes.