

Featured product!

Water removal issues?



Spiral Plug SP-V series "for water"

Fit!



almost **Double Speed**

Almost double* the evaporation rate of your aqueous samples

*When compared to the regular Spiral Plug (PTFE)



Recommended for those who need to evaporate **water!**

<Example of evaporation speed (5mL water evaporation time to dryness)>

•SP-P Plug (PTFE): 79minutes **SP-V Plug (Fluoroelastomer): 44minutes**

•Conditions : container: 20mL vial / set temperature: 40°C / exhaust speed 27L/min, 59%RH

Spiral Plug SP-V Series

Model No. SP-V□ *Add the plug size in □
(see the right table for the size)

Plug material Fluoroelastomer (FKM)

Inner diameter of container

Size Specification

4~7mm	V1
7~11mm	V2
11~17mm	V3
15~24mm	V4
24~32mm	V5

No.01

■Distributor

■Manufacturer

BC BioChromato

BioChromato, Inc.

1-12-19 Honcho, Fujisawa-City, Kanagawa, 251-0053 JAPAN

TEL: +81 466 238 382

Fax: +81 466 238 279

<https://www.biochromato.com/>

E-mail: oversea@bicr.co.jp