

Upgrade path to four position for your Smart Evaporator C1


4Kit for Smart Evaporator C1

“4kit” transforms your C1 into FOUR position simply by exchanging the head parts and heat block.



+

4 Kit



Smart Evaporator C 1

Smart Evaporator C1 with 4kit attached


4kit is a set of:

- 1 head part that accommodates four plugs (to be connected directly to the C1 plug adaptor)
- 4 each of P1 and P2 Spiral Plugs
- 4 sealing plugs
- 1 heat block that replaces the aluminum beads (one block can accommodate any of below vials)
- 4 each of 2mL, 4mL and 6mL vials



Four plugs



Heat block



Please check the type of your C1 in order to choose the fitting 4kit.

| Smart Evaporator C1 | 4Kit | |
|---|-----------------------|---------------------|
|  Type of attachable Smart Evaporator: Old Model (2012.3~2018.3) | Part number | 4kit UW4-C1 |
| | Heat block size | φ 50mm |
| | Compatible containers | 2mL ,4mL, 6mL vials |
| | Compatible samples | 4channels |
|  Type of attachable Smart Evaporator: New Model (2018.4~) | Part number | 4kit UW4-C1B |
| | Heat block size | φ 69mm |
| | Compatible containers | 2mL ,4mL, 6mL vials |
| | Compatible samples | 4channels |

■ Distributor

■ Manufacturer

BC BioChromato

BioChromato, Inc.

1-12-19Honcho, Fujisawa-City, Kanagawa, 251-0053 JAPAN
 TEL: +81-466-23-8382 Fax: +81-466-23-8279
<https://www.biochromato.com/> E-mail: oversea@bicr.co.jp

No.02